


# Carbonhand® glove size


6 sizes – XS, S, M, L, XL, XLW. Sizes slightly smaller than SEM Glove – fits smaller hands.



# Carbonhand<sup>®</sup> cord

## How to measure


- Cord in 3 lengths – 80, 90, 100 cm
- Measure from the top of the belt to the outer part of the shoulder bone (acromium) and from acromium to the wrist. Deduct 15 cm.
- If in doubt of which length to choose, take the longer one.



# Arm straps


If in doubt about 2 sizes, take the smaller one

## Upper Arm strap



Measured at the widest point

## Lower arm strap



Measured 140mm from the wrist line (5.5 inches)